

Antibodies Product Datasheet

Product Name:	KDM5D Antibody
Catalog Number:	C44383
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	KDM5D HY HYA JARID1D KIAA0234 SMCY
Protein Name:	Lysine-specific demethylase 5D (EC 1.14.11.-) (Histocompatibility Y antigen) (H-Y) (Histone demethylase JARID1D) (Jumonji/ARID domain-containing protein 1D) (Protein SmcY)
Other Name:	
Human Swiss Prot No.:	Q9BY66
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KDM5D Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year