

Antibodies Product Datasheet

Product Name:	KGP1 Antibody
Catalog Number:	C44385
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	PRKG1 PRKG1B PRKGR1A PRKGR1B
Protein Name:	cGMP-dependent protein kinase 1 (cGK 1) (cGK1) (EC 2.7.11.12) (cGMP-dependent protein kinase I) (cGKI)
Other Name:	
Human Swiss Prot No.:	Q13976
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KGP1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year