

Antibodies Product Datasheet

Product Name:	KHDR2 Antibody
Catalog Number:	C44386
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KHDR2
Concentration:	1 mg/ml
Gene Name:	KHDRBS2 SLM1
Protein Name:	KH domain-containing, RNA-binding, signal transduction-associated protein 2 (Sam68-like mammalian protein 1) (SLM-1) (hSLM-1)
Other Name:	
Human Swiss Prot No.:	Q5VWX1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KHDR2
Storage and Stability:	-20°C/1 year