

Antibodies Product Datasheet

Product Name:	KHDR3 Antibody
Catalog Number:	C44387
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	KHDRBS3 SALP SLM2
Protein Name:	KH domain-containing, RNA-binding, signal transduction-associated protein 3 (RNA-binding protein T-Star) (Sam68-like mammalian protein 2) (SLM-2) (Sam68-like phosphotyrosine protein)
Other Name:	
Human Swiss Prot No.:	O75525
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KHDR3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year