

Antibodies Product Datasheet

Product Name:	KI2L4 Antibody
Catalog Number:	C44394
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 280-360
Concentration:	1 mg/ml
Gene Name:	KIR2DL4 CD158D KIR103AS
Protein Name:	Killer cell immunoglobulin-like receptor 2DL4 (CD158 antigen-like family member D) (G9P) (Killer cell inhibitory receptor 103AS) (KIR-103AS) (MHC class I NK cell receptor KIR103AS) (CD antigen CD158d)
Other Name:	
Human Swiss Prot No.:	Q99706
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KI2L4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year