

Antibodies Product Datasheet

Product Name:	KI3X1 Antibody
Catalog Number:	C44398
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KI3X1
Concentration:	1 mg/ml
Gene Name:	KIR3DX1 KIR3DL0 LENG12
Protein Name:	KI3X1
Other Name:	Putative killer cell immunoglobulin-like receptor-like protein KIR3DX1 (Leukocyte receptor cluster member 12)
Human Swiss Prot No.:	Q9H7L2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.20% sodium azide.
Specificity:	This antibody detects endogenous levels of KI3X1 at Human
Storage and Stability:	-20°C/1 year