

Antibodies Product Datasheet

Product Name:	KIF5C Antibody
Catalog Number:	C44405
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIF5C
Concentration:	1 mg/ml
Gene Name:	KIF5C KIAA0531 NKHC2
Protein Name:	KIF5C
Other Name:	Kinesin heavy chain isoform 5C (Kinesin heavy chain neuron-specific 2)
Human Swiss Prot No.:	O60282
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.247% sodium azide.
Specificity:	This antibody detects endogenous levels of KIF5C at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year