

Antibodies Product Datasheet

Product Name:	KINH Antibody
Catalog Number:	C44407
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KINH
Concentration:	1 mg/ml
Gene Name:	KIF5B KNS KNS1
Protein Name:	KINH
Other Name:	Kinesin-1 heavy chain (Conventional kinesin heavy chain) (Ubiquitous kinesin heavy chain) (UKHC)
Human Swiss Prot No.:	P33176
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.266% sodium azide.
Specificity:	This antibody detects endogenous levels of KINH at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year