

Antibodies Product Datasheet

Product Name:	KKCC1 Antibody
Catalog Number:	C44413
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 400-480
Concentration:	1 mg/ml
Gene Name:	CAMKK1 CAMKKA
Protein Name:	Calcium/calmodulin-dependent protein kinase kinase 1 (CaM-KK 1) (CaM-kinase kinase 1) (CaMKK 1) (EC 2.7.11.17) (CaM-kinase IV kinase) (Calcium/calmodulin-dependent protein kinase kinase alpha) (CaM-KK alpha) (CaM-kinase kinase alpha) (CaMKK alpha)
Other Name:	
Human Swiss Prot No.:	Q8N5S9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KKCC1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year