

Antibodies Product Datasheet

Product Name:	KKLC1 Antibody
Catalog Number:	C44414
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	KKLC1 CXorf61
Protein Name:	Kita-kyushu lung cancer antigen 1 (KK-LC-1) (Cancer/testis antigen 83) (CT83)
Other Name:	
Human Swiss Prot No.:	Q5H943
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KKLC1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year