

Antibodies Product Datasheet

Product Name:	KLD8B Antibody
Catalog Number:	C44416
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLD8B
Concentration:	1 mg/ml
Gene Name:	KLHDC8B FP17659
Protein Name:	KLD8B
Other Name:	Kelch domain-containing protein 8B
Human Swiss Prot No.:	Q8IXV7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.56% sodium azide.
Specificity:	This antibody detects endogenous levels of KLD8B at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year