

Antibodies Product Datasheet

Product Name:	KLF16 Antibody
Catalog Number:	C44420
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	KLF16 BTEB4 NSLP2
Protein Name:	Krueppel-like factor 16 (Basic transcription element-binding protein 4) (BTE-binding protein 4) (Novel Sp1-like zinc finger transcription factor 2) (Transcription factor BTEB4) (Transcription factor NSLP2)
Other Name:	
Human Swiss Prot No.:	Q9BXX1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KLF16 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year