

Antibodies Product Datasheet

Product Name:	KLF9 Antibody
Catalog Number:	C44423
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLF9
Concentration:	1 mg/ml
Gene Name:	KLF9 BTEB BTEB1
Protein Name:	Krueppel-like factor 9 (Basic transcription element-binding protein 1) (BTE-binding protein 1) (GC-box-binding protein 1) (Transcription factor BTEB1)
Other Name:	
Human Swiss Prot No.:	Q13886
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KLF9
Storage and Stability:	-20°C/1 year