

Antibodies Product Datasheet

Product Name:	KLK5 Antibody
Catalog Number:	C44429
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	KLK5 SCTE UNQ570/PRO1132
Protein Name:	Kallikrein-5 (EC 3.4.21.-) (Kallikrein-like protein 2) (KLK-L2) (Stratum corneum tryptic enzyme)
Other Name:	
Human Swiss Prot No.:	Q9Y337
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KLK5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year