

Antibodies Product Datasheet

Product Name:	KLP6 Antibody
Catalog Number:	C44432
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLP6
Concentration:	1 mg/ml
Gene Name:	KLP6
Protein Name:	KLP6
Other Name:	Kinesin-like protein KLP6 (Kinesin-like protein 6)
Human Swiss Prot No.:	B7ZC32
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.275% sodium azide.
Specificity:	This antibody detects endogenous levels of KLP6 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year