

Antibodies Product Datasheet

Product Name:	KPCD3 Antibody
Catalog Number:	C44435
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	PRKD3 EPK2 PRKCN
Protein Name:	Serine/threonine-protein kinase D3 (EC 2.7.11.13) (Protein kinase C nu type) (Protein kinase EPK2) (nPKC-nu)
Other Name:	
Human Swiss Prot No.:	O94806
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KPCD3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year