

Antibodies Product Datasheet

Product Name:	KPYR Antibody
Catalog Number:	C44440
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	PKLR PK1 PKL
Protein Name:	Pyruvate kinase isozymes R/L (EC 2.7.1.40) (Pyruvate kinase 1) (R-type/L-type pyruvate kinase) (Red cell/liver pyruvate kinase)
Other Name:	
Human Swiss Prot No.:	P30613
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KPYR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year