

Antibodies Product Datasheet

Product Name:	KRA31 Antibody
Catalog Number:	C44449
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA31
Concentration:	1 mg/ml
Gene Name:	KRTAP3-1 KAP3.1 KRTAP3.1
Protein Name:	Keratin-associated protein 3-1 (High sulfur keratin-associated protein 3.1) (Keratin-associated protein 3.1)
Other Name:	
Human Swiss Prot No.:	Q9BYR8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KRA31
Storage and Stability:	-20°C/1 year