

Antibodies Product Datasheet

Product Name:	KRA41 Antibody
Catalog Number:	C44450
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA41
Concentration:	1 mg/ml
Gene Name:	KRTAP4-1 KAP4.10 KRTAP4-10 KRTAP4.1 KRTAP4.10
Protein Name:	Keratin-associated protein 4-1 (Keratin-associated protein 4-10) (Keratin-associated protein 4.1) (Keratin-associated protein 4.10) (Ultrahigh sulfur keratin-associated protein 4.10)
Other Name:	
Human Swiss Prot No.:	Q9BYQ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KRA41
Storage and Stability:	-20°C/1 year