

Antibodies Product Datasheet

Product Name:	KRA56 Antibody
Catalog Number:	C44453
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA56
Concentration:	1 mg/ml
Gene Name:	KRTAP5-6 KAP5.6 KRTAP5.6
Protein Name:	Keratin-associated protein 5-6 (Keratin-associated protein 5.6) (Ultrahigh sulfur keratin-associated protein 5.6)
Other Name:	
Human Swiss Prot No.:	Q6L8G9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KRA56
Storage and Stability:	-20°C/1 year