

Antibodies Product Datasheet

Product Name:	KRA59 Antibody
Catalog Number:	C44454
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA59
Concentration:	1 mg/ml
Gene Name:	KRTAP5-9 KAP5.9 KRN1 KRTAP5.9 UHSK1
Protein Name:	KRA59
Other Name:	
Human Swiss Prot No.:	P26371
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KRA59 at Human
Storage and Stability:	-20°C/1 year