

Antibodies Product Datasheet

Product Name:	KRA81 Antibody
Catalog Number:	C44455
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA81
Concentration:	1 mg/ml
Gene Name:	KRTAP8-1 KAP8.1 KRTAP8.1
Protein Name:	KRA81
Other Name:	
Human Swiss Prot No.:	Q8IUC2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KRA81 at Human/Mouse
Storage and Stability:	-20°C/1 year