

Antibodies Product Datasheet

Product Name:	KRT85 Antibody
Catalog Number:	C44468
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	KRT85 KRTHB5
Protein Name:	Keratin, type II cuticular Hb5 (Hair keratin K2.12) (Keratin-85) (K85) (Type II hair keratin Hb5) (Type-II keratin Kb25)
Other Name:	
Human Swiss Prot No.:	P78386
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KRT85 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year