

Antibodies Product Datasheet

Product Name:	KT33B Antibody
Catalog Number:	C44471
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	KRT33B HHA3-II HKA3B KRTHA3B
Protein Name:	Keratin, type I cuticular Ha3-II (Hair keratin, type I Ha3-II) (Keratin-33B) (K33B)
Other Name:	
Human Swiss Prot No.:	Q14525
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	KT33B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year