

Antibodies Product Datasheet

Product Name:	KY Antibody
Catalog Number:	C44479
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KY
Concentration:	1 mg/ml
Gene Name:	KY
Protein Name:	KY
Other Name:	Kyphoscoliosis peptidase (EC 3.4.-.-)
Human Swiss Prot No.:	Q8NBH2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.88% sodium azide.
Specificity:	This antibody detects endogenous levels of KY at Human/Mouse
Storage and Stability:	-20°C/1 year