

Antibodies Product Datasheet

Product Name:	LAP2B Antibody
Catalog Number:	C44486
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	TMPO LAP2
Protein Name:	Lamina-associated polypeptide 2, isoforms beta/gamma (Thymopoietin, isoforms beta/gamma) (TP beta/gamma) (Thymopoietin-related peptide isoforms beta/gamma) (TPRP isoforms beta/gamma) [Cleaved into: Thymopoietin (TP) (Splenin); Thymopentin (TP5)]
Other Name:	
Human Swiss Prot No.:	P42167
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LAP2B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year