

Antibodies Product Datasheet

Product Name:	LAPM5 Antibody
Catalog Number:	C44487
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LAPM5
Concentration:	1 mg/ml
Gene Name:	LAPTM5 KIAA0085
Protein Name:	Lysosomal-associated transmembrane protein 5 (Lysosomal-associated multitransmembrane protein 5) (Retinoic acid-inducible E3 protein)
Other Name:	
Human Swiss Prot No.:	Q13571
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LAPM5
Storage and Stability:	-20°C/1 year