

Antibodies Product Datasheet

Product Name:	LAR1B Antibody
Catalog Number:	C44488
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	LARP1B LARP2
Protein Name:	La-related protein 1B (La ribonucleoprotein domain family member 1B) (La ribonucleoprotein domain family member 2) (La-related protein 2)
Other Name:	
Human Swiss Prot No.:	Q659C4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LAR1B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year