

Antibodies Product Datasheet

Product Name:	LAT1L Antibody
Catalog Number:	C44491
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LAT1L
Concentration:	1 mg/ml
Gene Name:	SLC7A5P2 IMAA MMAA
Protein Name:	LAT1L
Other Name:	
Human Swiss Prot No.:	Q9GIP4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LAT1L at Human
Storage and Stability:	-20°C/1 year