

Antibodies Product Datasheet

Product Name:	LAT4 Antibody
Catalog Number:	C44493
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LAT4
Concentration:	1 mg/ml
Gene Name:	SLC43A2 LAT4 PP7664
Protein Name:	LAT4
Other Name:	Large neutral amino acids transporter small subunit 4 (L-type amino acid transporter 4) (Solute carrier family 43 member 2)
Human Swiss Prot No.:	Q8N370
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.167% sodium azide.
Specificity:	This antibody detects endogenous levels of LAT4 at Human/Mouse
Storage and Stability:	-20°C/1 year