

Antibodies Product Datasheet

Product Name:	LCE4A Antibody
Catalog Number:	C44501
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LCE4A
Concentration:	1 mg/ml
Gene Name:	LCE4A LEP8 SPRL4A
Protein Name:	LCE4A
Other Name:	
Human Swiss Prot No.:	Q5TA78
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LCE4A at Human
Storage and Stability:	-20°C/1 year