

Antibodies Product Datasheet

Product Name:	LCTL Antibody
Catalog Number:	C44503
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LCTL
Concentration:	1 mg/ml
Gene Name:	LCTL KLPB UNQ3022/PRO9820
Protein Name:	LCTL
Other Name:	Lactase-like protein (Klotho/lactase-phlorizin hydrolase-related protein)
Human Swiss Prot No.:	Q6UWM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.143% sodium azide.
Specificity:	This antibody detects endogenous levels of LCTL at Human/Mouse
Storage and Stability:	-20°C/1 year