

Antibodies Product Datasheet

Product Name:	LFTY1 Antibody
Catalog Number:	C44514
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LEFTY1 LEFTB LEFTYB UNQ278/PRO317
Protein Name:	Left-right determination factor 1 (Left-right determination factor B) (Protein lefty-1) (Protein lefty-B)
Other Name:	
Human Swiss Prot No.:	O75610
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LFTY1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year