

Antibodies Product Datasheet

Product Name:	LGAT1 Antibody
Catalog Number:	C44516
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 180-260
Concentration:	1 mg/ml
Gene Name:	LPGAT1 FAM34A KIAA0205
Protein Name:	Acyl-CoA:lysophosphatidylglycerol acyltransferase 1 (EC 2.3.1.-)
Other Name:	
Human Swiss Prot No.:	Q92604
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LGAT1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year