

Antibodies Product Datasheet

Product Name:	LHX6 Antibody
Catalog Number:	C44522
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LHX6
Concentration:	1 mg/ml
Gene Name:	LHX6 LHX6.1
Protein Name:	LIM/homeobox protein Lhx6 (LIM homeobox protein 6) (LIM/homeobox protein Lhx6.1)
Other Name:	
Human Swiss Prot No.:	Q9UPM6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LHX6
Storage and Stability:	-20°C/1 year