

Antibodies Product Datasheet

Product Name:	LIGO1 Antibody
Catalog Number:	C44524
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	LINGO1 LERN1 LRRN6A UNQ201/PRO227
Protein Name:	Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1 (Leucine-rich repeat and immunoglobulin domain-containing protein 1) (Leucine-rich repeat neuronal protein 1) (Leucine-rich repeat neuronal protein 6A)
Other Name:	
Human Swiss Prot No.:	Q96FE5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LIGO1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year