

Antibodies Product Datasheet

Product Name:	LIGO2 Antibody
Catalog Number:	C44525
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LINGO2 LERN3 LRRN6C UNQ9234/PRO31993
Protein Name:	Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 2 (Leucine-rich repeat neuronal protein 3) (Leucine-rich repeat neuronal protein 6C)
Other Name:	
Human Swiss Prot No.:	Q7L985
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LIGO2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year