

Antibodies Product Datasheet

Product Name:	LIPA1 Antibody
Catalog Number:	C44530
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 980-1060
Concentration:	1 mg/ml
Gene Name:	PPFIA1 LIP1
Protein Name:	Liprin-alpha-1 (LAR-interacting protein 1) (LIP-1) (Protein tyrosine phosphatase receptor type f polypeptide-interacting protein alpha-1) (PTPRF-interacting protein alpha-1)
Other Name:	
Human Swiss Prot No.:	Q13136
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LIPA1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year