

Antibodies Product Datasheet

Product Name:	LMF1 Antibody
Catalog Number:	C44538
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LMF1
Concentration:	1 mg/ml
Gene Name:	LMF1 C16orf26 TMEM112 HMFN1876 JFP11
Protein Name:	LMF1
Other Name:	Lipase maturation factor 1 (Transmembrane protein 112)
Human Swiss Prot No.:	Q96S06
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.144% sodium azide.
Specificity:	This antibody detects endogenous levels of LMF1 at Human/Mouse
Storage and Stability:	-20°C/1 year