

Antibodies Product Datasheet

Product Name:	LMTK1 Antibody
Catalog Number:	C44540
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1110-1190
Concentration:	1 mg/ml
Gene Name:	AATK AATYK KIAA0641 LMR1 LMTK1
Protein Name:	Serine/threonine-protein kinase LMTK1 (EC 2.7.11.1) (Apoptosis-associated tyrosine kinase) (AATYK) (Brain apoptosis-associated tyrosine kinase) (CDK5-binding protein) (Lemur tyrosine kinase 1) (p35-binding protein) (p35BP)
Other Name:	
Human Swiss Prot No.:	Q6ZMQ8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LMTK1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year