

Antibodies Product Datasheet

Product Name:	LNX1 Antibody
Catalog Number:	C44542
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LNX1
Concentration:	1 mg/ml
Gene Name:	LNX1 LNX PDZRN2 UNQ574/PRO1136
Protein Name:	E3 ubiquitin-protein ligase LNX (EC 6.3.2.-) (Ligand of Numb-protein X 1) (Numb-binding protein 1) (PDZ domain-containing RING finger protein 2)
Other Name:	
Human Swiss Prot No.:	Q8TBB1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LNX1
Storage and Stability:	-20°C/1 year