

Antibodies Product Datasheet

Product Name:	LONF1 Antibody
Catalog Number:	C44543
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LONF1
Concentration:	1 mg/ml
Gene Name:	LONRF1 RNF191
Protein Name:	LONF1
Other Name:	LON peptidase N-terminal domain and RING finger protein 1 (RING finger protein 191)
Human Swiss Prot No.:	Q17RB8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.350% sodium azide.
Specificity:	This antibody detects endogenous levels of LONF1 at Human
Storage and Stability:	-20°C/1 year