

Antibodies Product Datasheet

Product Name:	LONM Antibody
Catalog Number:	C44544
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LONM
Concentration:	1 mg/ml
Gene Name:	LONP1 PRSS15
Protein Name:	LONM
Other Name:	Lon protease homolog, mitochondrial (EC 3.4.21.-) (LONHs) (Lon protease-like protein) (LONP) (Mitochondrial ATP-dependent protease Lon) (Serine protease 15)
Human Swiss Prot No.:	P36776
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.251% sodium azide.
Specificity:	This antibody detects endogenous levels of LONM at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year