

Antibodies Product Datasheet

Product Name:	LPAR4 Antibody
Catalog Number:	C44547
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	LPAR4 GPR23 LPA4 P2RY9
Protein Name:	Lysophosphatidic acid receptor 4 (LPA receptor 4) (LPA-4) (G-protein coupled receptor 23) (P2Y purinoceptor 9) (P2Y9) (P2Y5-like receptor) (Purinergic receptor 9)
Other Name:	
Human Swiss Prot No.:	Q99677
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LPAR4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year