

Antibodies Product Datasheet

Product Name:	LPAR6 Antibody
Catalog Number:	C44548
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 70-150
Concentration:	1 mg/ml
Gene Name:	LPAR6 P2RY5
Protein Name:	Lysophosphatidic acid receptor 6 (LPA receptor 6) (LPA-6) (Oleoyl-L-alpha-lysophosphatidic acid receptor) (P2Y purinoceptor 5) (P2Y5) (Purinergic receptor 5) (RB intron encoded G-protein coupled receptor)
Other Name:	
Human Swiss Prot No.:	P43657
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LPAR6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year