

Antibodies Product Datasheet

Product Name:	LPCT4 Antibody
Catalog Number:	C44549
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 300-380
Concentration:	1 mg/ml
Gene Name:	LPCAT4 AGPAT7 AYTL3 LPEAT2
Protein Name:	Lysophospholipid acyltransferase LPCAT4 (EC 2.3.1.-) (1-acylglycerol-3-phosphate O-acyltransferase 7) (1-AGP acyltransferase 7) (1-AGPAT 7) (Acyltransferase-like 3) (Lysophosphatidylcholine acyltransferase 4) (Lysophosphatidylethanolamine acyltransferase
Other Name:	
Human Swiss Prot No.:	Q643R3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LPCT4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year