

Antibodies Product Datasheet

Product Name:	LPHN3 Antibody
Catalog Number:	C44550
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LPHN3 KIAA0768 LEC3
Protein Name:	Latrophilin-3 (Calcium-independent alpha-latrotoxin receptor 3) (CIRL-3) (Lectomedin-3)
Other Name:	
Human Swiss Prot No.:	Q9HAR2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LPHN3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year