

Antibodies Product Datasheet

Product Name:	LPP60 Antibody
Catalog Number:	C44553
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LPP60
Concentration:	1 mg/ml
Gene Name:	ASPG C14orf76
Protein Name:	60 kDa lysophospholipase (EC 3.1.1.5) [Includes: L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase); Platelet-activating factor acetylhydrolase (PAF acetylhydrolase) (EC 3.1.1.47)]
Other Name:	
Human Swiss Prot No.:	Q86U10
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LPP60
Storage and Stability:	-20°C/1 year