

Antibodies Product Datasheet

Product Name:	LPPR4 Antibody
Catalog Number:	C44555
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LPPR4 KIAA0455 PRG1
Protein Name:	Lipid phosphate phosphatase-related protein type 4 (EC 3.1.3.4) (Brain-specific phosphatidic acid phosphatase-like protein 1) (Plasticity-related gene 1 protein) (PRG-1)
Other Name:	
Human Swiss Prot No.:	Q7Z2D5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LPPR4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year