

Antibodies Product Datasheet

Product Name:	LRBA Antibody
Catalog Number:	C44556
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1970-2050
Concentration:	1 mg/ml
Gene Name:	LRBA BGL CDC4L LBA
Protein Name:	Lipopolysaccharide-responsive and beige-like anchor protein (Beige-like protein) (CDC4-like protein)
Other Name:	
Human Swiss Prot No.:	P50851
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LRBA Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year